

# Spectrometer Design Guide

A technical note

## Introduction


This guide provides some simple and easy to use design guidelines and formulas for designing, evaluating and comparing various diode array, diffraction grating based spectrometers designs. The input to the design process is the wavelength range you want to cover and the optical resolution by which you need to resolve the various structures in your spectrum (often peaks).

### CAUTION !

Spectrometer designs made by using this guide **should only be used as a starting point** in your design process. If you are going to implement a spectrometer in hardware you should always use a numerical simulation tool (for instance geometrical ray tracing) to make the final design.


## How a diode array spectrometer works

Basically, a spectrometer is an optical system consisting of two lenses/mirrors that produces an image of the input slit on the detector. In between the lenses/mirrors is placed a diffraction grating which disperses different wavelengths in different angles. This causes different wavelengths of light entering the input slit to be imaged to different position on the detector array.



On the following pages are shown two common spectrometer geometries; the transmission grating based and the crossed Czerny-Turner. Also, the figures defines the key design parameters of a spectrometer.

## Czerny-Turner




Minimum wavelength:	$\lambda_1$
Maximum wavelength:	$\lambda_2$
Wavelength range:	$\lambda_2 - \lambda_1$
Resolution:	$\Delta\lambda$
Center wavelength:	$\lambda_c = (\lambda_2 + \lambda_1)/2$

Angle of incidence:	$\alpha$
Diffraction angle:	$\beta$
	$\Phi = \beta - \alpha$

Grating groove density:	$G$
Focal length collimation:	$L_C$
Focal length focus:	$L_F$
Detector width:	$L_D$
Input slit width:	$\omega_{\text{slit}}$

## Transmission grating based




Minimum wavelength:	$\lambda_1$
Maximum wavelength:	$\lambda_2$
Wavelength range:	$\lambda_2 - \lambda_1$
Resolution:	$\Delta\lambda$
Center wavelength:	$\lambda_c = (\lambda_2 + \lambda_1)/2$

Angle of incidence:	$\alpha$
Diffraction angle:	$\beta$
	$\Phi = \beta - \alpha$

Grating groove density:	$G$
Focal length collimation:	$L_C$
Focal length focus:	$L_F$
Detector width:	$L_D$
Input slit width:	$\omega_{\text{slit}}$

## Flow chart



## The Spectrometer guide 8 steps

### Step 1: Choose geometry

The first step is to choose among the Czerny-Turner or LGL type geometries. For the Czerny Turner a typical value for  $\Phi$  is around  $30^\circ$  whereas transmission gratings are generally used in the Littrow configuration and 1<sup>st</sup> order where  $\alpha = \beta \Rightarrow \Phi = 0^\circ$ .

### Step 2: Choose grating

The second step is to choose a diffraction grating. Most grating vendors have an on-line catalogue where you can find one or more grating options to try in your design. You should choose a grating that has high diffraction efficiency in your wavelength. The important parameter that you shall use for the design in the next steps is the groove density  $G$ .

### Step 3: Calculate diffraction angle

The angle of incidence  $a$  on the grating and the diffraction angle  $b$  for the center wavelength  $\lambda_c$  are key parameters in the spectrometer design. These angles can be calculated once the grating groove density  $G$  and the total deflection  $\Phi$  is chosen.

### Step 4: Choose detector

The purpose of the spectrometer design is to disperse the wavelength range across the width of the detector array  $L_D$ . There are a large range of diode array detectors specifically designed for spectrometers. In general, if you need a compact spectrometer you should aim for a short detector (typically 1/4" or 6.4 mm). However, if you require a broad spectral range and/or a high resolution you should aim for a wide detector (typically 1/1" or 25.8 mm).

### Step 5: Calculate focal length of focus lens

Once the width of the detector is known you can calculate the focal length of the focusing mirror/lens.

### Step 6: Choose magnification

As mentioned earlier, the spectrometer is imaging the input slit to the detector and we generally want to have the slit as wide as possible to collect as much light through the input slit as possible. Therefore, the magnification in the system  $M$  should preferably be close to 1 which means that the width of the input slit ideally is imaged 1:1 onto the detector array.

### Step 7: Calculate focal length of collim lens

As in any imaging system the magnification is determined by the ratio between the focal lengths of the two lenses in the system. For a spectrometer this ratio has to be slightly modified due to the deflection along the beam path in the grating. However, once the magnification is chosen the focal length of the collimation mirror/lens can easily be calculated.

### Step 8: Calculate input slit width

The input slit width  $w_{\text{slit}}$  is determined by the required optical resolution  $\Delta\lambda$  and the magnification. Once you know your input slit width you are ready to evaluate if your spectrometer design is viable or you have to go back and change some of your choices for grating, detector, or magnification for instance.

## Evaluation of design

Once you have done a design iteration using the 8 steps described in the previous pages you should check whether this design is practical at all. Three things that you can easily check are the input slit width, diffraction limit of optics, and diffraction limit of grating as described below.

### Is the input slit width practical?

Input slits comes in widths down to 5 microns but such narrow slits will only allow very limited amount of light to enter your spectrometer. So, if your design requires a slit of 5 – 10 microns or less you could consider the following:

- loosening your requirement for the resolution
- choosing a wider detector and choosing a grating with a higher groove density

### Diffraction limit of optics

The formulas on the previous pages does not take into account that the optics can never produce a spot smaller than the diffraction limit. You can use the following formula to calculate the FWHM in wavelength of the smallest possible spot your optics can produce:

$$\Delta\lambda_{\text{diffraction}} = 1.028 \frac{\lambda_c M(\lambda_2 - \lambda_1)}{2 L_D \tan(\theta_{\text{NA}})}$$

If this value is larger than your required  $\Delta\lambda$ , your system is diffraction limited by the optics and you will not be able to obtain a better resolution than  $\Delta\lambda_{\text{diffraction}}$

### Diffraction limit of grating

The grating itself does also have a diffraction limited spot size (referred to as resolving power of the grating). The more grating lines that are being illuminated in a grating the better the resolution of the grating. The following formula gives the FWHM in wavelength of the smallest possible spot your grating can produce:

$$\Delta\lambda_{\text{diffraction}} = 0.84 \frac{\lambda_c \cos(\alpha)}{2 G L_G \tan(\theta_{\text{NA}})}$$

If this value is larger than your required  $\Delta\lambda$ , your system is diffraction limited by the optics and you will not be able to obtain a better resolution than  $\Delta\lambda_{\text{diffraction}}$

## About Ibsen Photonics

Ibsen was founded in 1991 by Per Ibsen under the name of Ibsen Micro Structures A/S. Today, 88% of Ibsen Photonics is majority owned by Foss A/S, a world leader in analytical solutions for the Food and Agricultural industries. Ibsen management and employees own 12 % of the shares in the company.

The Ibsen spirit combines the dynamic, entrepreneurial culture of a medium size company with a disciplined, operational mentality of a large corporation. With an average employee tenure of more than 10 years, Ibsen makes for a very effective organization that builds on more than 30 years of experience as a company.

Ibsen employs more than 75 people at our R&D and manufacturing facility in Denmark and achieved a turnover of more than 150 MDKK in 2023.

## Working with Ibsen Photonics

The core expertise of Ibsen Photonics lies in opto-mechanical design, grating technology and metrology. We master the cycle from optics, grating simulation and design, through optical and semiconductor production technologies, to high volume assembly, packaging and testing. Over the years we have developed many new designs, technologies and processes - many patented.

Our customers are large to medium-sized manufacturers of advanced optical devices and instruments, into which our products are integrated. With a highly organized production process, we are able to help customers obtain smooth instrument production, low unit-to-unit variation, high level of right first time, no field returns, and a low level of rework.

Our grating production facilities are world-class, including class 10 cleanroom facilities that we designed and built in 2000/2001, in which all environmental parameters are under continuous surveillance.

Our spectrometers are produced under strict quality control in our assembly facility in Denmark, certified to ISO 9001, ISO 13485, ISO 14001 and ISO 45001. This confirms Ibsen's capability to consistently produce high quality products that meet market standards and all regulatory requirements.

## Contact us

Ibsen Photonics A/S  
Ryttermarken 17  
DK-3520 Farum  
Denmark

**Telephone**  
**+45 4434 7000**



**Email**  
**[sales@ibsen.com](mailto:sales@ibsen.com)**

